




Instruction Manual for Ion Drop Counter




Product No.: ZC1036

■ Description

This product consists of an Ion Drop Counter, syringe set, accessories and sensor electrode. It is used to explore the acid-base drop experiment and determinate the concentration of some ionized solution.

■ Product Specification

| No. | Name | Picture | Quantity | Remark |
|-----|------------------|---|----------|---|
| 1 | Ion Drop Counter |  | 1 set | |
| 2 | Syringe set |  | 1 set | 1 set of syringe, 2 sets of three-way valve, 1 drip nozzle |
| 3 | Accessories |  | 1 set | 1 piece of rotary union, 1 piece of aluminum rod, and 1 piece of binder clip |

| | | | | |
|---|------------------|---|----------|--|
| 4 | Plum screw |  | 1 piece | |
| 5 | Hand Screw |  | 3 pieces | |
| 6 | Sensor electrode |  | 1 set | |

Tab.1

■ **Product General Assembly Figure**





Fig.1: Ion Drop Counter

(The iron stand and the temperature sensor probe are provided separately)

■ **Product General Assembly Procedures**

Assembly is not required, as shown in Fig. 1.

■ Function Area



■ Specification of Sensor Electrode

| Channel | Name | Function | Measuring Range | Resolution | Accuracy | Other specification/ first and second value for calibration |
|---------|----------------------------------|----------------------------|--|------------|-------------|---|
| 1 | Drop counter | Count the number of drops | 0~∞c/0~∞mL | 1c/0.1mL | ±1c/ ±0.1mL | Maximum drop rate: 30 drops/s |
| 2 | Universal input port | Connect to universal input | Refer to the specification of connected sensor | | | |
| 3 | pH sensor | | 0 ~14 | 0.01 | ±0.2pH | 4.00:9.18 |
| | Nitrite ion sensor | | 0.5ppm~ 4600ppm | 0.4ppm | | 46ppm: 4600ppm (0.1mol/L:0.001mol/L) |
| | Nitrate ion sensor | | 0.6ppm~ 6200ppm | 0.4ppm | | 62ppm: 6200ppm (0.1mol/L:0.001mol/L) |
| | Chloride ion sensor | | 1.8ppm~ 3550ppm | 0.3ppm | | 35ppm: 3550ppm (0.1mol/L:0.001mol/L) |
| | Ammonium ion sensor | | 0.9ppm~ 1800ppm | 0.1ppm | | 18ppm: 1800ppm (0.1mol/L:0.001mol/L) |
| | Sodium ion sensor | | 0.2ppm~ 2300ppm | 0.2ppm | | 23ppm: 2300ppm (0.1mol/L:0.001mol/L) |
| | Potassium ion sensor | | 0.2ppm~ 39000ppm | 3ppm | | 390ppm: 39000ppm (0.1mol/L:1mol/L) |
| | Calcium ion sensor | | 0.4ppm~ 4000ppm | 0.7ppm | | 40ppm: 4000ppm (0.1mol/L:0.001mol/L) |
| | Dissolved CO ₂ sensor | | 0.2ppm~ 440ppm | 0.1ppm | | 44ppm: 440ppm (0.1mol/L:0.001mol/L) |
| 4 | Temperature sensor | | -40°C~+135°C | 0.1°C | ±0.6°C | |

Note: The electrode in the standard configuration includes only pH sensor electrode. Other electrodes are provided separately.

■ Typical Experiment

Acid-base neutralization drop

Heat by acid-base reaction

Produce nitrogen dioxide by decomposition of nitrite

Detect concentration of sodium ion, potassium and calcium ions

■ Notes

1. Before experiment, connect sensors to sensor electrode.
2. The experiment should be carried out by using the specific experiment software.
3. The other ion sensor electrode are provided separately. You can purchase them if required.

■ Maintenance

The device should be placed in a cardboard box if not in use. Use a damp cloth to wipe out the dust on the surface lightly. Avoid scratches caused by sharp articles on the surface.